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
Taking the underlying concepts described for the desktop simulation one step further is the recreation of the real physical patient in a realistic physical clinical environment. This is done because real patients don't live in a virtual world – they and the clinicians who care for them work in an actual physical world with real people. Thus, computerized mannequin stands in for the patient, and a variety of equipment can be used (either real clinical equipment or computer-driven replicas) to monitor and treat the patient.


The mannequin-based simulator has a computer representation of the patient similar to that in a desktop simulator, replacing the videos, drawings, and animations with actual functions of the “plastic person.” Among the functions that these mannequin-based simulators can replicate are:

- Simulation Modalities - Mannequin-based Simulation
- Spontaneous breathing (and the ability to breathe for the patient with a bag or ventilator)
- Real-time display of electronically monitored information (e.g. ECG, oxygen saturation, etc.)
- Pulses, heart sounds, breath sounds, pupil size, pupil response to light
- Obstruction of various parts of the airway

Mannequin-based Patient Simulation

พจนานุกรม : http://cisl.stanford.edu/what_is/sim_modalities/mannequin_sim.html

CAE Fidelis™ Lucina childbirth simulator

All the stages of labor and delivery
10 simulated clinical experiences (SCEs)

- A normal delivery
- An instrumental vaginal delivery
- Fetal tachycardia due to maternal pyrexia
- Shoulder dystocia
- Major post-partum hemorrhage due to uterine atony
- Breech delivery
- Fetal central nervous system depression by narcotics given to the mother
- Major post-partum hemorrhage
- Maternal cardio-respiratory arrest
- Eclampsia
- Umbilical cord prolapse

Instructor's Workstation with Muse

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


FETUS

4. Power on delivery fetus and lubricate.
5. Attach umbilical cord to delivery fetus.
6. Install delivery fetus.
Presentation: Vertex
Loading Position: Occiput posterior (face up).

ABDOMEN

7. Connect umbilical cord to placenta.
8. Insert intact placenta in placenta pouch. Attach pouch to pelvis.
9. Install delivery abdomen (connect speakers and palpable contraction system).

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1. Install rotation ring.
2. Use standard birth canal (lubricate).
3. Install uterine funnel with delivery cervix (lubricate).
4. Power on delivery fetus and lubricate.
5. Attach umbilical cord to delivery fetus.
6. Install delivery fetus. Presentation: Vertex Loading Position: Occiput posterior (face up).
7. Connect umbilical cord to placenta.
8. Insert intact placenta in placenta pouch. Attach pouch to pelvis.
9. Install delivery abdomen (connect speakers and palpable contraction system).

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Synopsis

You have been called to a residence of a 65-year-old patient who complains he has had increasing shortness of breath, fatigue and weight gain over the past couple of weeks. The dyspnea has suddenly worsened overnight. He experienced difficulty when tying his shoes this morning and states his pants are snug around the waist.

This SCE has five states, which are automatically or manually progress as described below.

During the initial assessment in **State 1 Initial Presentation**, the patient presents with increasing dyspnea, HR in the 110s, BP in the 220s/140s, RR in the 40s and SpO2 in the mid 90s on room air. Auscultation of breath sounds reveals crackles in all fields. The patient remains conscious and alert. The learner is expected to administer oxygen and initiate cardiac monitoring and intravenous therapy. If the learner applies oxygen the SCE automatically advance to **State 2 Slight Improvement with Oxygen**.

In **State 2 Slight Improvement with Oxygen**, HR is in the 90s, BP in the 190s/110s, RR in the 30s and SpO2 in the mid 90s on 2 LPM oxygen via nasal cannula. If the learner fails to administer oxygen within four minutes, end the SCE. If the learner unsuccessfully manages the patient, it is recommended to repeat the simulation until a positive outcome occurs. If Nitrates are administered, the SCE is manually transitioned to **State 3 Slight Improvement with Nitrates**.

In **State 3 Slight Improvement with Nitrates**, the patient manifests increasing dyspnea. The patient's HR remains in the 120s, BP is 120s/90s, RR in the upper 30s and SpO2 drops to the high 80s. The learner is expected to obtain a 12-lead ECG, administer sublingual nitrates and consider ventilatory assistance. In 120 seconds the SCE automatically advances to **State 4 Continued Improvement**.

In **State 4 Continued Improvement**, the patient's HR slows to the 70s, BP is 140s/90s, RR is in the mid teens and SpO2 improves to the mid 90s. The learner is expected to administer additional nitrates as needed, consider ACE Inhibitors and modify oxygen delivery as indicated. If appropriate therapy continues, in 120 seconds the SCE automatically advances to **State 5 Transport to Facility**.

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Norma L. Female
33 years old
70.0 kg
Postpartum

Content Management

- Overview
- Background**
- Preparation
- Notes

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Operating Mode | Vaginal Delivery

- Prepartum and Latent Phase
- Active Phase
- Vaginal Delivery
- Cesarean Section
- Postpartum
- Non-gravid Female

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Cervix | Closed

- Closed
- One Finger (0%)
- One Finger (50%)
- Two Fingers (70%)
- 2cm (70%)
- 3cm (90%)
- 4cm (90%)
- 5cm (90%)

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Postpartum Initial Uterine State | Contracted

- Contracted
- Boggy
- Fully Inverted
- Partially Inverted

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Norma L. Female
33 years old
70.0 kg
Postpartum

Content Management

Overview Background
Preparation Notes

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Save
Save As
Save A Copy

Rename Scenario

New State
Rename Selected State
Delete Selected State

ALL CONDITIONS
Cardiovascular
Hemorrhage
Neurologic
Respiratory
Interventions

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
Sounds
Button

Obstetrical Basic

- Rate of Descent: 1.0
- Fetal State: Active Sleep
- Contraction Frequency: 4
- Contraction Amplitude: 60
- Contraction Duration: 80
- Patient Pushing: Off
- Early Deceleration Magnitude: None
- Late Deceleration Amplification Factor: 1.0
- Shoulder Dystocia: Disabled
- Shoulder Dystocia Resolution: None
- Extraction of Posterior Arm:
- Arrested Labor: Disable
- Arrested Labor Trigger Station: 0
- Arrested Labor Resolved by Traction: No

Heart Basic

- Blood Pressure: Modeled
- Heart Rate: Modeled
- Heart Rate Factor: 1.00
- Cardiac Rhythm: Modeled
- Arterial Catheter: Peripheral Artery
- Central Venous Catheter: Right Atrium
- PA Catheter: Pulmonary Artery
- PA Balloon: Deflated
- Defib:
- Pacing Current: 0
- Pacing Rate: 80
- Pacing Capture Threshold: 50
- Cold Fluid Inject:

Head

Apply To Both Eyes: Off

- 5 mm
- Auto
- Eyes: Blink Speed: Normal
- Reactive pupils: Yes
- Light Reactivity Speed: Brisk
- Seizures: Off
- ICP: Modeled
- NMB: Modeled
- Temperature: Body: 36.5
- Temperature: Blood: 37.0

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ALL CONDITIONS

- Cardiovascular
- Hemorrhage
- Neurologic
- Respiratory
- Interventions

| Cardiovascular | Respiratory | Hemorrhage | Neurologic |
|----------------------------|---------------------|-----------------------------|-------------------|
| BP: Hypertension | Acidosis | Bleeding Channels Off/On | ICP Elevation |
| BP: Hypotension | Alkalosis | Bleeding Location Channel 1 | Temperature Blood |
| Cardiac Tamponade | Apnea | Bleeding Location Channel 2 | Temperature Body |
| Heart Rate: Bradycardia | Bradypnea | Bleeding Type Channel 1 | |
| Heart Rate: Tachycardia | Chest Wall Dynamics | Bleeding Type Channel 2 | |
| Patient Stability | CO2: Hypercapnea | Hemorrhage: Blood | |
| Ventricular Failure: Both | CO2: Hypocapnea | Hemorrhage: Plasma | |
| Ventricular Failure: Left | Desaturation | | |
| Ventricular Failure: Right | Intrapleural Volume | | |
| | I to E Ratio | | |
| | Tachypnea | | |

Lung Basic

- Bronchial Occlusion: Off
- Respiratory Rate: Modeled
- Respiratory Rate Factor: 1.00
- Shunt Fraction: 0.02
- SPO2: Modeled
- NMB: Modeled
- Tidal Volume: Modeled
- Intrapleural Vol: Left: 0
- Intrapleural Vol: Right: 0
- Fraction of Inspired O2: 21

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Essential Of simulation

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Writing a Simulation

Please fill in the table PRIOR to working with the simulator:

State = Write in the Name that is obvious to anyone who may be using this e.g., Primary Assessment, Patient Feels Breat hless

A = Fill in any Airway related features that need to be programmed e.g. Swollen tongue

B = Fill in any Breathing related features that need to be programmed e.g., Respiratory rate, SpO₂, breath sounds, etc.

C = Fill in any Circulation related features that need to be programmed e.g., Heart rate, blood pressure, capillary refill, urinary output, etc.

D = Fill in any Disability (Neurological) related features that need to be programmed e.g., Eyes – blink speed, seizures, increased intracranial pressure

E = Fill in any Exam related features that need to be programmed e.g., Temperature (remember program in body & blood – with Blood 0.5°C higher), bowel sounds, trauma features not covered by B/C

Events = Things that simulator currently cannot do e.g., Movement of extremities

Lab results/Transitions = If you are using them and what they are

Example Template

| State | A | B | C | D | E | Events | Transitions |
|-------|---|---|---|---|---|--------|-------------|
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |

(Wilford, 2011)

Worked Example

Theme: Asthma for newly qualified/graduated nurses, assessment and initial intervention

Aim of session: For newly qualified/graduated nurses to assess, recognize and treat a patient with asthma in an acute medical unit

Learning Outcome:

- To demonstrate how to assess a patient in respiratory distress using ABCDE Framework
- To use communication skills to reassure and explain to the patient as required
- To administer nebulizer as per local protocol and assess the effectiveness of treatment

Name of file: Asthma Recognition

Patient History:

John Devine is a 40-year-old man admitted yesterday evening with exacerbation of asthma due to a chest infection. He usually takes a Salbutamol Inhaler 2 puffs 4 times a day. He does not smoke. Upon admission, the patient had a temperature of 38.2°C, RR 28, HR 105 irregular, BP 100/70, SpO₂ 91% on room air. Blood cultures and a sputum specimen were taken and he was started on prophylactic antibiotics, steroids and nebulizers (see drug chart). It is now 7 am the next morning. Please assess and treat John as necessary.

Props/Equipment:

Simulator, drug chart with appropriate drugs including the prescribed nebulizers, medical and nursing notes, observation/assessment chart, fluid chart, selection of oxygen masks and salbutamol inhaler. Patient to be semi-prone on 4 LPM oxygen via nasal cannula

Debrief points:

ABCDE

Prompt treatment – Discuss what if nebulizer did not work. Also, discuss Chronic Obstructive Pulmonary Disease (COPD)

Communication – Explore

Give reference list for current protocol of asthma

(Wilford, 2011)

Worked Example

| State | A | B | C | D | E | Events | Transitions |
|------------------------------------|-----|--|---|------------------------|---------------------------------|--|---|
| First meeting at 7am | nil | RR- 25-28 SpO ₂ on 4 LPM - 92% Bilateral wheeze | HR 110 irreg - ectopics BP 140/90 Caprefill > 3 secs | Eyes - Blink Fast | T Body = 37.5 T Blood = 38.0 | Anxious Sweating, pale Blood sugar = 9 mmols/L or 162 mg/dL Finding it hard to talk | nil |
| Deteriorates if no treatment given | nil | RR- 28-33 SpO ₂ on 4 LPM - 89% Bilateral wheeze - Louder | HR 120 irreg - ectopics BP 100/70 Caprefill > 3 secs | Eyes - Blink Fast | T Body = 37.5 T Blood = 38.0 | More anxious Sweating, pale Blood sugar = 9 mmols/L or 162 mg/dL | nil |
| Nebulizer given | nil | RR- 22 over 3 minutes SpO ₂ increases to 95% over 3 minutes Bilateral wheeze - Less | HR 130 (effects of Beta 2) BP 110/80 Caprefill < 3 secs | Eyes - Blink Fast | T Body = 37.5 T Blood = 38.0 | Anxious Sweating, less pale Blood sugar = 9 mmols/L or 162 mg/dL | After 3 minutes, move to recovery automatically |
| Recovery | nil | RR- 20 over 1 minute SpO ₂ - 95% Now wheeze | HR 100 with ectopics BP 120/80 Caprefill > 3 secs | Eyes - Blink Normal | T Body = 37.5 T Blood = 38.0 | Anxious Sweating, perfused skin Blood sugar = 9 mmols/L or 162 mg/dL | nil |

(Wilford, 2011)

Appendix B - SCE Development

Writing a simulation

Theme: _____

Aim of session:

Learning Outcome:

(Include at least 1 technical and 1 non-technical, no more than 3-5)

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


1. Prepartum and Latent Phase
2. Active Phase
3. Vaginal Delivery
4. Cesarean Section
5. Postpartum
6. Non-gravid Female

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Conditions Guidelines for Programming MFS Lucina

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Respiratory: Desaturation

Respiratory: Respiratory Rate

| Tachypnea | | Bradypnea | |
|--------------|----|--------------|----|
| Reset | 11 | Reset | 11 |
| Increased | 15 | Increased | 10 |
| Elevated | 18 | Intermediate | 9 |
| Borderline | 20 | Mild | 7 |
| Intermediate | 22 | Moderate | 6 |
| Mild | 25 | Severe | 5 |
| Moderate | 28 | Profound | 3 |
| Severe | 31 | Extreme | 2 |
| Profound | 33 | | |
| Extreme | 36 | | |

Respiratory: Desaturation

| Desaturation | SpO ₂ Value |
|--------------|------------------------|
| Reset | 98% |
| High 90s | 96-97% |
| Mid 90s | 94-96% |
| Low 90s | 91-93% |
| High 80s | 87-90% |
| Mid 80s | 84-86% |
| Low 80s | 80-83% |
| High 70s | 77-80% |
| Mid 70s | 74-77% |
| Low 70s | 69-71% |
| Less than 70 | <69% |

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Cardiovascular: Blood Pressure

| Hypertension | | Hypotension | |
|----------------|-----------|----------------|----------|
| Reset | 110s/70s | Reset | 110s/70s |
| Increased | 120s/80s | Decreased | 100s/70s |
| Pre-Borderline | 130s/80s | Pre-Borderline | 100s/60s |
| Borderline | 140s/90s | Borderline | 90s/50s |
| Mild | 150s/90s | Mild | 80s/40s |
| Moderate | 160s/100s | Moderate | 70s/40s |
| Severe | 170s/100s | Severe | 60s/30s |
| Profound | 190s/110s | Profound | 50s/30s |
| Extreme | 220s/120s | Extreme | 40s/30s |

Cardiovascular: Heart Rate

| Tachycardia | | Bradycardia | |
|----------------|-----------|----------------|----------|
| Reset | 70s | Reset | 70s |
| Increased | High 70s | Decreased | Mid 60s |
| Elevated | 80s | Pre-Borderline | Low 60s |
| Pre-Borderline | 90s | Borderline | Mid 50s |
| Borderline | 100s | Intermediate | Low 50s |
| Intermediate | 110s | Mild | High 40s |
| Mild | 120s | Moderate | Mid 40s |
| Moderate | 130s | Severe | Low 40s |
| Severe | 140s | Extreme | Mid 30s |
| Supra | 150s | Acute | Low 30s |
| Profound | 160s | | |
| Extreme | 170s | | |
| Acute | High 170s | | |

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

| Ventilation Mode or Simulator Feature | Works with Lucina | Operation / Adjustments | Limitations / Expected Performance |
|--|-------------------|---|---|
| Controlled Mechanical Ventilation (CMV) | Yes | Adjustments made on the ventilator | V_T max: 500 mL No official support for PEEP |
| Ventilator Triggering (i.e. SIMV mode) | Yes | Adjustments made on the ventilator | Requires flow or pressure trigger at minimum setting |
| Self-Refilling Ventilation Bag (i.e. Bag-valve-mask (BVM)) | Yes | | Mask leak during assisted/controlled ventilation |
| Non-Refilling Reservoir Bag (i.e. Mapleson circuits) | Limited | Controlled ventilation | Minimal spontaneous exhalation – therefore Lucina will NOT fully self-refill reservoir bag |
| Reservoir/Breathing Bag in Circle System (i.e. Anesthesia machine) | Limited | Controlled ventilation | Minimal spontaneous exhalation – therefore Lucina will NOT fully self-refill reservoir bag |
| Oxygen Delivery (FiO_2) (via ventilator, face mask, nasal cannula) | Yes | Requires SW setting of FiO_2 for physiologic response | Lucina does NOT automatically respond to ventilator FiO_2 settings |
| Ventilator Applied Positive End-Expiratory Pressure (PEEP) | Limited* | Requires SW setting of PEEP for physiologic response | *Limited by ventilator V_T Maximum PEEP at set V_T in CMV mode before over-distention: V_T ~400 mL: PEEP 5 cmH ₂ O |
| Exhaled CO ₂ | No | | |
| Capnogram | Limited | TouchPro ONLY | Lucina will NOT demonstrate capnogram on a clinical monitor |
| Volatile Anesthetic Agent Functionality (i.e. Isoflurane) | Limited | SW (Virtual) agent ONLY | Lucina does NOT detect anesthetic agents |
| Variable Bronchial Resistance | No | SW control of the hardware | Bronchial Occlusion ONLY (L/R, both) - CAN demonstrate increased peak inspiratory pressure (PIP) |
| Variable Lung Compliance | No | SW setting ONLY | CANNOT demonstrate increased plateau pressure (and PIP) while on the ventilator |
| Difficult Airway Capabilities | No | | Does NOT have swollen tongue, airway occluder or laryngospasm - CANNOT perform cricothyrotomy |

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

Initial Assessment

Parameters Notes Custom Variable

| | |
|--------------------------------------|-----------------------|
| Contraction Frequency | 5 ctx/10 min over 240 |
| Contraction Amplitude | 80 mmHg over 240 sec. |
| Temperature: Blood | 37.5 °C |
| Temperature: Body | 37 °C |
| Heart Rate Factor | 1.15 |
| Systemic Vascular Resistance Factor | 1.35 |
| Venous Capacity Factor | 0.9 |
| Systemic Arteries Compliance Factor | 0.9 |
| Pulmonary Arteries Compliance Factor | 0.9 |
| Respiratory Rate Factor | 1.75 |
| Shunt Fraction | 0.08 |
| Tidal Volume Factor | 0.9 |

Transitions

Create

IF Dilatation > 9 cm
GOTO Four Hours Later
ELSE

The patient is a 27-year-old primigravida who was admitted in labor four hours ago, at 40 weeks gestation. Mother's blood group is O positive. Group B streptococcus at 36 weeks was negative. State 2 is four hours later, however, the delivery is programmed to occur over 20 minutes in real time.

This SCE consists of three states that are transitioned manually at the facilitator's discretion, except for State 1 which automatically transitions to State 2 Four Hours Later when fetus is at +1 station (approximately four minutes).

State 1 Initial Assessment

The patient demonstrates a HR in the 80s, BP in the 110s/70s, RR in the mid 20s and SpO₂ in the upper 90s on room air and a temperature of 37.0C.

Vaginal examination on admission revealed that the cervix is 100% effaced and 7 cm dilated.

The fetus is in cephalic presentation (vertex) at -2, and the membranes were found to be ruptured upon vaginal examination.

The fundal height is 35 cm and fetal heart rate (FHR) is 140s beats per minute (bpm). The patient is presenting regular contractions at four-minute intervals. No analgesia has been requested.

The learner is expected to provide support and information to patient and family, recognize cervical dilatation, evaluate FHR and contractions and recognize regular contraction frequency.

CAE Healthcare : MFS Lucina Normal Delivery (SCE)

Four Hours Later

Parameters Notes Custom Variable

| | |
|-------------------------------------|---------------------------|
| Contraction Frequency | 6 ctx/10 min over 120 sec |
| Contraction Amplitude | 100 mmHg over 120 sec. |
| Heart Rate Factor | 1.43 over 30 sec. |
| Systemic Vascular Resistance Factor | 1.5 over 30 sec. |
| Venous Capacity Factor | 0.75 over 30 sec. |
| Venous Return Resistance Factor | 0.8 over 30 sec. |
| Respiratory Rate Factor | 2 over 15 sec. |

Transitions

State 2 Four Hours Later

The patient demonstrates a HR in the 80s, BP in the 120s/80s, RR in the upper 20s, SpO2 in the upper 90s on room air and the FHR is in the 140s.

The patient states that she feels the urge to push.

Facilitator should activate the Patient Pushing Parameter.

The learner is expected to complete a vaginal examination and find the cervix 10 cm dilated and fetus at +1 station, recognize the start of second stage labor, ensure the presence of appropriate personnel in the room, position and prepare the patient for delivery, and complete the delivery of the fetus. Facilitator should manually transition to State 3 after delivery of the fetus.

Delivery of Placenta

Parameters Notes Custom Variable

| | |
|-------------------------------------|-------------------|
| Heart Rate Factor | 1.05 over 30 sec. |
| Systemic Vascular Resistance Factor | 1.35 over 30 sec. |
| Venous Capacity Factor | 0.9 over 30 sec. |
| Respiratory Rate Factor | 1.2 over 15 sec. |

Transitions

State 3 Delivery of Placenta


The patient demonstrates a HR in the 70s, BP in the 110s/70s, RR in the teens and SpO2 in the upper 90s on room air.

The learner is expected to deliver and examine the placenta.

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The patient is a 22-year-old primiparous woman who was admitted in labor eight hours ago, at 39 weeks gestation.

There is no prior history of medical complications, and her pregnancy has been uneventful with routine and regular prenatal care. Mother's blood group is O positive.

Group B streptococcus is negative.

Epidural analgesia and continuous external FHR monitoring have been in place for four hours.

This SCE consists of eleven states that automatically transition or are transitioned manually at the facilitator's discretion, based on the learner's interventions.

If the facilitator chooses, the shoulder dystocia can be relieved by one of the internal and/or external maneuvers and the SCE will automatically advance to State 8 Dystocia Resolved.

If the facilitator chooses, the internal and/or external maneuvers can be unsuccessful in resolving the shoulder dystocia and the facilitator can advance the scenario to State 10 Zavanelli Maneuver, if performed by the learner, and the SCE then automatically advances to State 11 Emergency Cesarean Section when the Zavanelli maneuver is performed.

**CAE Healthcare : MFS Lucina
Shoulder Dystocia (SCE)**

Initial Assessment

Parameters Notes Custom Variable

| | |
|--------------------------------------|--------------|
| Shoulder Dystocia | Enabled |
| Contraction Amplitude | 81 mmHg |
| Contraction Frequency | 5 ctx/10 min |
| Temperature: Body | 37 °C |
| Temperature: Blood | 37.5 °C |
| Heart Rate Factor | 0.93 |
| Venous Capacity Factor | 1.05 |
| Systemic Arteries Compliance Factor | 0.55 |
| Pulmonary Arteries Compliance Factor | 0.5 |
| Respiratory Rate Factor | 0.9 |

Transitions Create

IF Shoulder Dystocia is Yes
GOTO Arrest of Descent Due to Shoulder Dystocia

Arrest of Descent Due to Shoulder Dystocia

Parameters Notes Custom Variable

| | |
|-----------------|--------------------|
| FHR Signal Loss | 100 % over 60 sec. |
|-----------------|--------------------|

Transitions Create

State 1 Initial Assessment

The patient demonstrates a HR in the 80s, BP in the 100s/50s, RR in the low 10s and SpO2 in the upper 90s on room air. Vaginal examination demonstrates that the cervix is 100% effaced and 10 cm dilated.

The fetus is in the cephalic presentation (vertex) at +5 and membranes are ruptured. The fundal height is 35 cm. The fetal head delivers, but the shoulders remain inside the maternal pelvis in left occiput transverse (LOT).

The fetal heart rate (FHR) is 140 beats per minute (bpm) with moderate variability, accelerations and decelerations. Uterine contractions are every two to three minutes, lasting 60-90 seconds and are firm to palpation.

The SCE automatically advances to State 2 Arrest of Descent Due to Shoulder Dystocia after delivery of the fetal head.

State 2 Arrest of Descent Due to Shoulder Dystocia

The learner is expected to recognize shoulder dystocia and complete the internal and/or external maneuvers for the resolution of shoulder dystocia.

If the facilitator chooses, shoulder dystocia is relieved by one of the internal and/or external maneuvers.

NOTE: The learner is expected to progress through the appropriate sequence of maneuvers - State 3 to State 7.

The facilitator selects which maneuver, if any, will resolve shoulder dystocia. The SCE automatically transitions to State 8 Dystocia Resolved when the selected resolution is performed.

McRoberts Maneuver Resolves Dystocia

Parameters Notes Custom Variable

| | |
|-------------------------------------|-------------------|
| Shoulder Dystocia Resolution | McRoberts |
| Heart Rate Factor | 1.02 over 30 sec. |
| Venous Capacity Factor | 1 over 30 sec. |
| Systemic Arteries Compliance Factor | 0.5 over 30 sec. |
| Respiratory Rate Factor | 1 over 15 sec. |

Transitions Create

IF Shoulder Dystocia is No

GOTO Dystocia Resolved

State 3 McRoberts Maneuver Resolves Dystocia

The patient demonstrates a HR in the 80s, BP in the 110s/50s, RR in the low teens, SpO2 in the 90s and FHR in the 140s with moderate variability.

The learner is expected to ask assistants to use the McRoberts maneuver to resolve the shoulder dystocia.

Suprapubic Pressure Resolves Dystocia

Parameters Notes Custom Variable

| | |
|-------------------------------------|---------------------|
| Shoulder Dystocia Resolution | Suprapubic Pressure |
| Heart Rate Factor | 1.02 over 30 sec. |
| Venous Capacity Factor | 1 over 30 sec. |
| Systemic Arteries Compliance Factor | 0.5 over 30 sec. |
| Respiratory Rate Factor | 1 over 15 sec. |

Transitions Create

IF Shoulder Dystocia is No

GOTO Dystocia Resolved

State 4 Suprapubic Pressure Resolves Dystocia

The patient demonstrates a HR in the 80s, BP in the 110s/50s, RR in the low teens, SpO2 in the 90s and FHR in the 140s with moderate variability.

The learner is expected to ask an assistant to apply suprapubic pressure to resolve the shoulder dystocia.

Rotation of the Posterior Shoulder Resolves...

| | | |
|-------------------------------------|--------------------------------|-------------------|
| Parameters | Notes | + Custom Variable |
| Shoulder Dystocia Resolution | Rotation of Posterior Shoulder | |
| FHR Baseline | 105 bpm over 60 sec. | |
| FHR Variability Coefficient | 0.65 over 30 sec. | |
| Heart Rate Factor | 1.02 over 30 sec. | |
| Venous Capacity Factor | 1 over 30 sec. | |
| Systemic Arteries Compliance Factor | 0.5 over 30 sec. | |
| Respiratory Rate Factor | 1 over 15 sec. | |

Transitions Create

IF Shoulder Dystocia is No

GOTO Dystocia Resolved

Rotation of the Anterior Shoulder Resolves ...

| | | |
|-------------------------------------|-------------------------------|-------------------|
| Parameters | Notes | + Custom Variable |
| Shoulder Dystocia Resolution | Rotation of Anterior Shoulder | |
| FHR Baseline | 105 bpm over 60 sec. | |
| FHR Variability Coefficient | 0.65 over 30 sec. | |
| Heart Rate Factor | 1.02 over 30 sec. | |
| Venous Capacity Factor | 1 over 30 sec. | |
| Systemic Arteries Compliance Factor | 0.5 over 30 sec. | |
| Respiratory Rate Factor | 1 over 15 sec. | |

Transitions Create

IF Shoulder Dystocia is No

GOTO Dystocia Resolved

State 5 Rotation of the Posterior Shoulder Resolves Dystocia

The patient demonstrates a HR in the 80s, BP in the 110s/50s, RR in the low teens, SpO2 in the 90s and FHR in the 100s with variable decelerations and minimal variability.

The learner is expected to rotate the posterior shoulder to resolve the shoulder dystocia.

State 6 Rotation of the Anterior Shoulder Resolves Dystocia

The patient demonstrates a HR in the 80s, BP in the 110s/50s, RR in the low teens, SpO2 in the 90s and FHR in the 100s with variable decelerations and minimal variability.

The learner is expected to rotate the anterior shoulder to resolve the shoulder dystocia.

Extraction of the Posterior Arm

Parameters Notes Custom Variable

| | |
|-------------------------------------|----------------------|
| FHR Baseline | 105 bpm over 60 sec. |
| FHR Variability Coefficient | 0.65 over 30 sec. |
| Heart Rate Factor | 1.02 over 30 sec. |
| Venous Capacity Factor | 1 over 30 sec. |
| Systemic Arteries Compliance Factor | 0.5 over 30 sec. |
| Respiratory Rate Factor | 1 over 15 sec. |

Transitions Create

IF Shoulder Dystocia is No
GOTO Dystocia Resolved

Dystocia Resolved

Parameters Notes Custom Variable

Transitions Create

Delivery of Placenta

Parameters Notes Custom Variable

| | |
|-------------------------------------|-------------------|
| Heart Rate Factor | 0.95 over 30 sec. |
| Systemic Vascular Resistance Factor | 1.1 over 30 sec. |
| Venous Capacity Factor | 0.9 over 30 sec. |
| Systemic Arteries Compliance Factor | 0.75 over 30 sec. |
| Respiratory Rate Factor | 1.2 over 30 sec. |

Transitions Create

State 7 Extraction of Posterior Arm

The patient demonstrates a HR in the 80s, BP in the 110s/50s, RR in the low teens, SpO2 in the 90s and FHR in the 100s with variable decelerations and minimal variability. The learner is expected to extract the posterior arm to resolve the shoulder dystocia and complete a vaginal delivery.

NOTE:

The facilitator should select 'Extract Now' under the parameter Extraction of Posterior Arm (Fetal and Labor Parameters, Basic parameters) when this maneuver is performed.

State 8 Dystocia Resolved

The vaginal delivery is completed.

State 9 Delivery of Placenta

The patient demonstrates a HR in the 70s, BP in the 110s/60s, RR in the mid teens and SpO2 in the upper 90s on room air. The learner is expected to deliver and examine the placenta.

**CAE Healthcare : MFS Lucina
Shoulder Dystocia (SCE)**

Zavanelli Maneuver

Parameters Notes Custom Variable

| | |
|-------------------------------------|---------------------|
| Shoulder Dystocia Resolution | Zavanelli |
| FHR Baseline | 95 bpm over 60 sec. |
| FHR Variability Coefficient | 0.35 over 30 sec. |
| Venous Capacity Factor | 1 over 30 sec. |
| Systemic Arteries Compliance Factor | 0.5 over 30 sec. |
| Respiratory Rate Factor | 1 over 15 sec. |

Transitions Create

IF Shoulder Dystocia is No
GOTO Emergency Cesarean Section

Emergency Cesarean Section

Parameters Notes Custom Variable

Transitions Create

State 10 Zavanelli Maneuver

performed, the patient demonstrates a HR in the 80s, BP in the 110s/50s, RR in the low teens, SpO2 in the 90s and FHR in the 90s with variable decelerations and absent variability.

The learner is expected to replace the fetal head and immediately transfer the patient to the operating room.

The SCE automatically transitions to State 11 Emergency Cesarean Section when the Zavanelli maneuver is performed.

State 11 Emergency Cesarean Section

The patient demonstrates a HR in the 80s, BP in the 110s/50s, RR in the low teens, SpO2 in the 90s and FHR is 90s with prolonged deceleration and absent variability.

The learner is expected to assemble the equipment required, prepare the patient for emergency cesarean section and complete the delivery via cesarean section.

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After Blood Products

Treatment with Prostaglandins

Treatment with Bimanual Uterine Compression

Treatment with Intrauterine Balloon

Uterine Artery Embolization

Surgery

State 1 Initial Assessment

, the patient demonstrates a HR in the 120s, BP in the 90s/60s, RR in the mid 20s and SpO2 in the upper 90s on room air. The patient starts to complain of dizziness and loss of vision.

The learner is expected to assess the state of uterine contractions, place a urinary catheter, start external uterine massage, start intravenous access draw labwork, initiate crystalloid infusions at a high rate, lower the patient's head, and initiate oxygen and continuous ECG monitoring.

| Initial Assessment | |
|--------------------------------------|-----------------|
| Parameters | Notes |
| Uterine Massage | Recurring atony |
| Fluid Loss Blood | 800 mL |
| Temperature: Body | 37 °C |
| Temperature: Blood | 37.5 °C |
| Heart Rate Factor | 1.3 |
| Systemic Vascular Resistance Factor | 1.15 |
| Venous Capacity Factor | 0.95 |
| Systemic Arteries Compliance Factor | 0.85 |
| Pulmonary Arteries Compliance Factor | 0.85 |
| Respiratory Rate Factor | 1.7 |
| Shunt Fraction | 0.07 |
| Functional Residual Capacity | 2000 mL |

Transitions

Uterus Contracts

Parameters Notes Custom Variable

| | |
|--------------------------------------|------------------|
| Uterine Massage | Effective |
| Vaginal Bleeding Severity | None |
| Heart Rate Factor | 1.2 over 30 sec. |
| Left Ventricle Contractility Factor | 1.1 over 30 sec. |
| Right Ventricle Contractility Factor | 1.1 over 30 sec. |
| Venous Capacity Factor | 0.8 over 30 sec. |
| Systemic Arteries Compliance Factor | 0.9 over 30 sec. |
| Pulmonary Arteries Compliance Factor | 0.9 over 30 sec. |
| Respiratory Rate Factor | 1.2 over 15 sec. |

Transitions

State 2 Uterus Contracts

the patient's condition demonstrates a HR in the 100s, BP in the 100s/70s, RR in the upper teens, and SpO2 in the upper 90s on room air.

The learner is expected to consider the need for colloids. If blood products are administered, the facilitator should manually transition to State 4 After Blood Products.

Boggy Uterus

Parameters Notes Custom Variable

| | |
|--------------------------------------|------------------|
| Uterine Massage | Recurring atony |
| Vaginal Bleeding Severity | Moderate |
| Systemic Vascular Resistance Factor | 1.3 over 30 sec. |
| Systemic Arteries Compliance Factor | 0.9 over 30 sec. |
| Pulmonary Arteries Compliance Factor | 0.9 over 30 sec. |
| Respiratory Rate Factor | 1.61 |

Transitions

State 3 Boggy Uterus

The patient's condition demonstrates a HR in the 130s, BP in the 80s/60s, RR in the mid 20s, and SpO2 in the upper 90s on room air.

The uterus is again in a boggy state. The learner is expected to reassess the state of uterine contraction, initiate external uterine massage, reassess the cause of bleeding and consider the need for colloids. If blood products administered, the facilitator should manually transition to State 4 After Blood Products.

After Blood Products

Parameters Notes Custom Variable

| | |
|--------------------------------------|-------------------|
| PRBC Infusion | 250 mL |
| Whole Blood Infusion | 250 mL |
| Colloid Infusion | 500 mL |
| Heart Rate Factor | 1.15 over 30 sec. |
| Left Ventricle Contractility Factor | 1.05 over 30 sec. |
| Right Ventricle Contractility Factor | 1.05 over 30 sec. |
| Respiratory Rate Factor | 1.2 |

Transitions

Treatment with Prostaglandins

Parameters Notes Custom Variable

| | |
|-------------------------------------|------------------|
| Vaginal Bleeding Severity | Mild |
| Systemic Vascular Resistance Factor | 1.2 over 30 sec. |
| Venous Capacity Factor | 0.9 over 30 sec. |

Transitions

Treatment with Bimanual Uterine Compression

Parameters Notes Custom Variable

| | |
|-------------------------------------|------------------|
| Vaginal Bleeding Severity | Mild |
| Systemic Vascular Resistance Factor | 1.2 over 30 sec. |
| Venous Capacity Factor | 0.9 over 30 sec. |

Transitions

State 4 After Blood Products

The patient's condition demonstrates a HR in the 80s-100s, BP in the 100s-110s/60s-70s, RR in the upper teens, SpO2 in the upper 90s.

The learner is expected to assess the patient's response to the administration of blood products.

State 5 Treatment with Prostaglandins

, the patient's condition demonstrates a HR in the 110s, BP in the 90s/60s, RR in the upper teens, SpO2 in the upper 90s. The uterus remains boggy.

The learner is expected to administer prostaglandins as ordered and evaluate their effect on uterine contraction.

State 6 Treatment with Bimanual Uterine Compression

The patient's condition demonstrates a HR in the 110s, BP in the 90s/60s, RR in the upper teens, SpO2 in the upper 90s. The uterus remains boggy.

The learner is expected to initiate bimanual uterine compression and evaluate the efficacy of interventions on uterine contraction.

Treatment with Intrauterine Balloon —

Parameters Notes + Custom Variable

| | |
|-------------------------------------|------------------|
| Vaginal Bleeding Severity | Mild |
| Systemic Vascular Resistance Factor | 1.2 over 30 sec. |
| Venous Capacity Factor | 0.9 over 30 sec. |

Transitions

State 7 Treatment with Intrauterine Balloon

The patient's condition demonstrates a HR in the 110s, BP in the 90s/60s, RR in the upper teens, SpO2 in the upper 90s.

The learner is expected to insert or assist with the insertion of an intrauterine balloon, fill balloon with 500 mL fluid, and evaluate the effectiveness of intervention on uterine contraction.

Uterine Artery Embolization —

Parameters Notes + Custom Variable

| | |
|--------------------------------------|----------------|
| Vaginal Bleeding Severity | None |
| Systemic Arteries Compliance Factor | 1 over 30 sec. |
| Pulmonary Arteries Compliance Factor | 1 over 30 sec. |

Transitions

State 8 Uterine Artery Embolization

The patient's condition demonstrates a HR in the 100s, BP in the 90s/60s-70s, RR in the upper teens, SpO2 in the upper 90s.

The learner is expected to prepare the patient for transfer to Interventional Radiology for embolization.

State 9 Surgery,

The patient's condition demonstrates a HR in the 50s, BP in the 80s/40s, RR less than 10s, SpO2 in the upper 80s.

The learner is expected to prepare the patient for emergency surgery.

Surgery —

Parameters Notes + Custom Variable

| | |
|--------------------------------------|-------------------|
| Fluid Loss Blood | 500 mL |
| Baroreceptor Gain (Cardiac) Factor | 0.75 over 45 sec. |
| Heart Rate Factor | 0.17 over 45 sec. |
| Systemic Vascular Resistance Factor | 1.15 over 45 sec. |
| Venous Capacity Factor | 0.85 over 45 sec. |
| Systemic Arteries Compliance Factor | 0.8 over 45 sec. |
| Pulmonary Arteries Compliance Factor | 0.8 over 45 sec. |
| Respiratory Rate Factor | 0.3 over 15 sec. |
| Shunt Fraction | 0.25 |

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Postpartum



Major Postpartum Hemorrhage Due to Uterine Atony

Make sure the mannequin is setup as follows:

BIRTHING MECHANISM

FETUS

ABDOMEN

1. Remove rotation ring.

4. Ensure no fetus is installed.

5. Install prepartum abdomen.

2. Connect and install the bleeding tank. Let the placenta pouch rest over it.

3. Install postpartum birth canal with boggy/contracted uterus.

Connect pneumatic and bleeding tubes.

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Ineffective - Massage of the uterus does not cause it to contract.

Effective - With **Effective**, successful compression of the uterus will fully contract the boggy uterus. Five effective compressions are required before contraction occurs. The feel of the uterus will not change for the first four massages, i.e., the uterus stays fully boggy for massages one to four. With the fifth massage, the uterus fully contracts over approximately 90 seconds and the message "Boggy uterus was declared resolved" appears in the event log.

A single massage/compression consists of pressing, holding, and releasing the uterus. With each successful massage, the message "Uterine compression was applied" is logged in the event log, independent of uterus mode (**Effective** or **Recurring atony**).

In order to completely contract the uterus, one must execute five successful compressions in **Effective** or **Recurring atony** mode.

Additionally, when **Vaginal Bleeding Severity** is active (Mild, Moderate, Severe, or Profuse), each massage will produce a gush of blood in addition to the bleeding enabled by the vaginal bleeding severity parameter (see **Vaginal Bleeding** in the *Using Fidelis Lucina* section of this guide). This gush of blood occurs every time one presses on the uterus, whether or not the message "Uterine compression was applied" is logged in the event log, and is independent of the uterus mode (**Effective** or **Recurring atony**).

Recurring atony - With **Recurring atony**, the function is similar to **Effective** (i.e., the boggy uterus contracts after five successful logged massages). However, 90 seconds after contraction the uterus returns to the boggy state, and the message "Boggy uterus was declared resolved" does not appear in the event log. In this mode, if massage is interrupted for 90 seconds or longer before it fully contracts, the massage count is reset, even when the uterine compressions are marked as effective (i.e., logged as "Uterine compression was applied").

For example, if the learner applies three logged massages then waits 90 seconds or longer to apply the next massage, the 'massage count' resets. From this point, five effective consecutive massages (i.e., without interruption) are required before the uterus contracts.

หมายเหตุ : it will take approximately 60 seconds for the uterus to reach the boggy state. After six successful external massages have been recorded in the Event Log, the message "Boggy Uterus was Declared Resolved" will appear in the Event Log. At this time the boggy uterus will begin deflating which takes approximately 90 seconds, and the learner will then be able to palpate the contracted uterus. The contracted uterus will begin to become boggy again in approximately 90 seconds.

K10