CAE Fidelis™ Lucina

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A Vaginal Breech Delivery
Peggy Williams


Age: 22 years old
Weight: 61.0 kg
Base: Peggy Williams
Operating Mode: Vaginal Delivery

Overview
The patient is a 22-year-old Gravida 2, Para 1 who was admitted in labor 15 minutes ago, at 38 weeks gestation. Mother's blood group is B positive. Group B streptococcus at 38 weeks was negative. Vaginal examination on admission revealed cervix 100% effaced and fully dilated. The fetus is in uncertain presentation at +3 and the membranes are ruptured and positive for ferning. Fundal height is 32 cm and the fetal heart rate (FHR) is in the 150s beats per minute (bpm). The patient is presenting with regular and painful contractions. No monitoring is in place yet, and the patient begins to push. This delivery is programmed to occur over 16 minutes in real time.

This SCE consists of three states that are transitioned manually at the facilitator's discretion.

In State 1 Initial Assessment, the patient demonstrates a HR in the 80s, BP in the 120s/70s, RR in the upper 20s and SpO2 in the upper 90s on room air, a temperature of 36.8C and FHR in the 150s. Vaginal discharge is clear and non-odorous. Contractions are occurring every five minutes, 30-seconds duration, of moderate intensity upon palpation. Fetal heart tones are auscultated and traced above umbilicus. A vaginal exam reveals 100% effaced and fully dilated, fetal skin folds palpable, unable to palpate vertex landmarks. Findings are discussed with the interprofessional team. The learner is expected to assess the patient's presenting condition, discuss health history, establish labor pattern and fetal presentation, noting the possible breech presentation. The learner is expected to analyze the data and information gathered in this state and communicate with the interdisciplinary team. Admission orders are required and fetal presentation is to be confirmed with ultrasound.

In State 2 Breech Delivery, the contraction pattern is every two minutes, 60 to 80 second duration and moderate to strong intensity and FHR in the 150s. Maternal vital signs are BP in the 120/70s, HR in the 80s, RR

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muse

Search

Major Postpartum Hemorrhage Due to U...
Donna Chapman

Maternal Cardiopulmonary Arrest
Leticia Coyle

Non-gravid Female
Norma L. Female

Pregnant Female
Stennette Montobe

Shoulder Dystocia
Julie Thornton

Umbilical Cord Prolapse
Emma Mason

Page 2 / 2

Open Library New SCE...

Shoulder Dystocia
Julie Thornton

Age: 22 years old
Weight: 86.0 kg
Base: Julie Thornton
Operating Mode: Vaginal Delivery

View as PDF
Add To Favorites

Overview
The patient is a 22-year-old primiparous woman who was admitted in labor eight hours ago, at 39 weeks gestation. There is no prior history of medical complications, and her pregnancy has been uneventful with routine and regular prenatal care. Mother's blood group is O positive. Group B streptococcus is negative. Epidural analgesia and continuous external fetal monitoring have been in place for four hours.

This SCE consists of eleven states that automatically transition or are transitioned manually at the facilitator's discretion, based on the learner's interventions. If the facilitator chooses, the shoulder dystocia can be relieved by one of the internal and/or external maneuvers and the SCE will automatically advance to **State 9 Dystocia Resolved**. If the facilitator chooses, the internal and/or external maneuvers can be unsuccessful in resolving the shoulder dystocia and the facilitator can advance the scenario to **State 10 Zavanelli Maneuver**, if performed by the learner, and the SCE then automatically advances to **State 11 Emergency Cesarean Section** when the Zavanelli maneuver is performed.

In **State 1 Initial Assessment**, the patient demonstrates a HR in the 80s, BP in the 100s/50s, RR in the low 10s and SpO2 in the upper 90s on room air. Vaginal examination demonstrates that the cervix is 100% effaced and 10 cm dilated. The fetus is in the cephalic presentation (vertex) at +3 and membranes are ruptured. The fundal height is 35 cm. The fetal head delivers, but the shoulders remain inside the maternal pelvis in left occiput transverse (LOT). The fetal heart rate (FHR) is 140 beats per minute (bpm) with moderate variability, accelerations and decelerations. Uterine contractions are every two to three minutes, lasting 60-90 seconds and are firm to palpation.

The SCE automatically advances to **State 3 Arrest of Descent Due to Shoulder Dystocia** after delivery of the

Review Run

muse admin

Mannequin Setup

Vaginal Delivery - Vertex

Vaginal Delivery

Make sure the mannequin is setup as follows:

BIRTHING MECHANISM

FETUS

ABDOMEN

1. Install rotation ring.
2. Use standard birth canal (lubricate).
3. Install uterine funnel with delivery cervix (lubricate).
4. Power on delivery fetus and lubricate.
5. Attach umbilical cord to delivery fetus.
6. Install delivery fetus.
Presentation: Vertex
Loading Position: Occiput posterior (face up).
7. Connect umbilical cord to placenta.
8. Insert intact placenta in placenta pouch.
Attach pouch to pelvis.
9. Install delivery abdomen (connect speakers and palpable contraction system).

Cancel Continue

The Mannequin Setup Screen - Vaginal Delivery Vertex Presentation

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The screenshot shows the mUSE software interface. On the left is a search bar and a list of simulation scenarios. The main window displays the title 'Major Postpartum Hemorrhage Due to Uterine Atony' for a patient named Donna Chapman. The patient's details include: Age 39 years old, Weight 63.0 kg, Base Donna Chapman, and Operating Mode Postpartum. An overview section describes the patient's history and the simulation's states. At the bottom, there are 'Review' and 'Run' buttons.

The diagram, titled 'The Mannequin Setup Window', illustrates the correct setup for a postpartum simulation. It shows three views: 'BIRTHING MECHANISM', 'FETUS', and 'ABDOMEN'. The 'FETUS' view shows a red 'X' over the fetus, indicating it should not be installed. The 'ABDOMEN' view shows the correct placement of the prepartum abdomen. A list of instructions is provided below the diagrams:

1. Remove rotation ring.
2. Connect and install the bleeding tank. Let the placenta pouch rest over it.
3. Install postpartum birth canal with boggy/contracted uterus. Connect pneumatic and bleeding tubes.
4. Ensure no fetus is installed.
5. Install prepartum abdomen.

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